Get the Blood Drop - Palm

ALTERNATIVE SITE TESTING (AST) PALM

Use the **clear** AST endcap.*

1. Insert the lancet as described in *Prepare the Lancing Device*.
2. Attach and lock the **clear** AST endcap. Twist the white band until you hear or feel a click and the white and black raised marks line up as shown.

“Click” to lock

CAUTION: Do not use the **clear** AST endcap on your fingertips.

3. Select a puncture site from a fleshy area on the palm. Avoid veins, moles, bones, and tendons.


4. Pull back the re-set handle to prepare the lancing device. The release button is raised when the lancing device is ready.
5. Press the lancing device firmly against the puncture site, then press the release button UNTIL THE DEVICE CLICKS.


(If the device does not puncture, reset by repeating Step 4.)

6. Do not take the device away from your palm. Maintain steady pressure until a small, round blood drop forms under the clear endcap.

7. Lift the lancing device straight up and away from the skin without smearing the blood drop.
8. Test immediately.

Do not test the blood sample from palm if you get:

- Smears blood
- Runny blood
- Clotted blood
- Clear fluid mixed with the blood

9. Wipe the puncture site with a clean, dry tissue and maintain pressure on the site until the bleeding stops.

See: **Eject / Dispose of the Lancet**

Cleaning and Disinfection

Clean and disinfect your lancing device once a week. **Use only Super Sani-Cloth germicidal wipes** containing 55% isopropyl alcohol and 0.5% ammonium chloride. The wipes are proven safe to use with the lancing device.

Cleaning involves the removal of visible dirt and debris, but does not reduce the risk for transmission of infectious diseases.

Disinfecting (if performed properly) reduces the risk of transmitting infectious diseases.

Warning: *If your lancing device is being operated by a second person who is providing testing assistance to you, the lancing device should be disinfected prior to use by the second person.*

Super Sani-Cloth germicidal wipes are available for purchase online. For more information, contact Customer Service at 1-800-268-7200.

Do not use bleach (sodium hypochlorite).

For more information:

Health Canada advisory. Blood lancing devices for personal use may transmit blood-borne viruses if used on more than one person. Government of Canada; 3/25/09. <http://www.healthycanadians.gc.ca/recall-alert-rappel-avis/hc-sc/2009/13372a-eng.php>


Signs of Deterioration: Your lancing device has been tested for 156 cycles of the cleaning and disinfection procedure (one cycle per week for 3 years). The lancing device has been demonstrated to withstand 3 years of cleaning and disinfection without damage. However, if your device malfunctions for any reason, contact Customer Service at 1-800-268-7200 to troubleshoot and provide a solution.

Cleaning Your Lancing Device

Supplies Needed for Cleaning:

- Super Sani-Cloth germicidal wipes

1. Clean the lancing device with Super Sani-Cloth germicidal wipes until all soil is removed.
2. Allow to air dry.


Disinfecting Your Lancing Device

Supplies Needed for Disinfecting:


- Super Sani-Cloth germicidal wipes
- Timing device

1. Before disinfecting, clean your lancing device as described in *Cleaning Your Lancing Device*.
2. Remove the endcap and eject the lancet, if present.
3. Get a new Super Sani-Cloth germicidal wipe.
4. **Set timer for 2 minutes and start timing.**
5. **LANCING DEVICE BODY:** Wipe the body continuously. **You must wipe all surfaces and keep them wet for 2 minutes.**

Continue to step 6

6. **Set timer for 2 minutes and start timing again.**
7. **ENDCAP:** Wipe the endcap continuously. **You must wipe the inside and outside and keep them wet for 2 minutes.**

8. Allow lancing device body and endcap to air dry.

Customer Service

Warranty: Your lancing device is warranted for 3 years.

If you need assistance or have questions about blood sugar testing, please contact Customer Service at 1-800-268-7200.

MICROLET NEXT lancing device made in Poland.

Manufactured for:
Ascensia Diabetes Care Holdings AG
Peter Merian-Strasse 90
4052 Basel, Switzerland
1-800-268-7200

TM see www.ascensiadibetes.ca/tm-mc

Ascensia Diabetes Care Holdings AG
Peter Merian-Strasse 90
4052 Basel, Switzerland

For information on patents and related licenses, see:
www.patents.ascensia.com

Ascensia, the Ascensia Diabetes Care logo, and Microlet are trademarks and/or registered trademarks of Ascensia Diabetes Care Holdings AG.


All other trademarks are properties of their respective owners and are used solely for informative purposes. No relationship or endorsement should be inferred or implied.

© 2018 Ascensia Diabetes Care Holdings AG. All rights reserved.

90002101 / 90002102 Rev. 07/18

Microlet®
next Lancing Device

INSTRUCTIONAL INSERT

- Lancet ejector for safe lancet disposal.
- **Black/white** adjustable endcap for fingertip testing only.
- **Clear** AST endcap for alternative site testing only.*
- The MICROLET®NEXT lancing device and MICROLET® lancets are used to draw capillary blood from the fingertip or palm for blood glucose testing.

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⚠ WARNING


- Ask your health care professional if Alternative Site Testing (AST) is right for you.
- Do not calibrate a continuous glucose monitoring device from an AST result.
- Do not calculate an insulin dose based on an AST result.

Do not use AST under the following conditions:

- If you think your blood sugar is low (hypoglycaemia).
- When blood sugar is changing rapidly (after a meal, insulin dose, or exercise).
- If you are unable to feel symptoms of low blood sugar (hypoglycaemic unawareness).
- If you get alternative site blood sugar results that do not agree with how you feel.
- During illness or times of stress.
- If you will be driving a car or operating machinery.

POTENTIAL BIOHAZARD

- The MICROLET®NEXT lancing device is intended for single patient use only. It should not be used for assisted blood draws by health care providers or at health care provision sites, and should never be shared with anyone else, even a family member, due to risk of infection.
- Do not reuse lancets. Used lancets are not sterile. Use a new lancet each time you test.
- Always dispose of used test strips and lancets as medical waste or as advised by your health care professional.
- Always wash your hands with soap and water and dry them well before and after testing or handling the meter, lancing device, or test strips.

* May not be included in all kits. Available on request from **Customer Service**.

Prepare the Lancing Device

Warning: Keep out of reach of children. This kit contains small parts that could cause suffocation if accidentally swallowed.

1. Wash hands and the puncture site with soap and warm water. Rinse and dry well.


2. Unlock the endcap. Twist the white band so the white raised mark and black indented mark line up as shown.


3. Pull off the endcap.


4. Loosen the round protective cap on the lancet by rotating the cap ¼ turn, but **do not remove the cap**.

5. Insert the lancet firmly into the lancing device. **Push** until the lancet comes to a complete stop.

6. Carefully twist off the round protective cap. Save the cap for disposing of the used lancet.

Get the Blood Drop - Fingertip

FINGERTIP TESTING

Use the **black/white** adjustable endcap. Do not use the clear AST endcap.*


1. Attach the **black/white** adjustable endcap.
2. Lock the endcap. Twist the white band until you hear or feel a click and the white and black raised marks line up as shown.


3. Rotate the endcap dial to adjust the puncture depth. The amount of pressure applied to the puncture site also affects puncture depth.


4. Pull back the re-set handle to prepare the lancing device. The release button is raised when the lancing device is primed.


5. Press the lancing device firmly against the **side** of the finger and press the release button UNTIL THE DEVICE CLICKS.

6. Stroke your hand and finger toward the puncture site to form a drop of blood. Do not squeeze around the puncture site.

7. Test immediately after you have formed a small, round blood drop, as shown.

8. Wipe the puncture site with a clean, dry tissue and **maintain pressure** on the site until the bleeding stops.


Eject / Dispose of the Lancet

Always dispose of the used lancet as medical waste to prevent injury or contamination to others. Use a new lancet each time you test. **Do not reuse lancets.**


Your MICROLET NEXT lancing device has a lancet ejector. Do not use your fingers to remove the lancet from the lancing device.


1. Carefully remove the endcap from the lancing device. Twist the white band of the endcap so the white raised mark and the black indented mark line up. Pull off the endcap.

2. Place the round protective cap that you saved onto a flat surface. With the used lancet still in the lancing device, push the lancet into the center of the cap.

3. Point the device toward a container for sharp objects (an empty plastic bottle with a lid, for example). Label the container properly and keep the container out of the reach of children.

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